

# CONTENTS

Title	page
ABOUT LAYERTEC / ABOUT THIS CATALOG	02
CONTINUOUS WAVE AND NANOSECOND LASER OPTICS	03
· Introduction to continuous wave and nanosecond laser optics	03
· Components for F <sub>2</sub> Lasers	04
· Components for ArF Lasers	06
· Components for KrF, XeCl and XeF Lasers	08
· Components for Ti:Sapphire lasers	10
· Components for diode lasers	12
· Broadband and scanning mirrors	14
· Components for Yb:YAG, Yb:KGW and Yb-doped fiber lasers	16
· High power Nd:YAG/Nd:YVO <sub>4</sub> laser optics	18
· Components for the harmonics of Nd:YAG and Nd:YVO <sub>4</sub> lasers	20
· Components for other Nd:YAG or Nd:YVO <sub>4</sub> wavelengths	22
· Components for Ho:YAG and Tm:YAG lasers	24
· Components for Er:YAG lasers and the 3µm region	26
· Components for optical parametrical oscillators (OPOs)	28
PICOSECOND LASER OPTICS	30
FEMTOSECOND LASER OPTICS	32
· Introduction to femtosecond laser optics	32
· Femtosecond laser optics for Ti:Sapphire and related lasers	34
– Standard femtosecond laser optics (bandwidth ~120nm)	34
– Broadband femtosecond laser optics (bandwidth ~300nm)	36
– Ultra broadband femtosecond laser optics (bandwidth ~1octave)	38
– Silver mirrors for femtosecond lasers	40
· Femtosecond laser optics optimized for third order dispersion	42
· High power femtosecond laser optics	44
· Overview about laser induced damage thresholds of fs laser optics	45
· Components for the second harmonic of the Ti:Sapphire laser	46
· Components for higher harmonics of the Ti:Sapphire laser	48
· Gires-Tournois-Interferometer (GTI) mirrors	50
· Optics for femtosecond lasers in the 1100 – 1600nm wavelength range	52

Title	page
FILTERS FOR LASER APPLICATIONS	54
LOW LOSS OPTICAL COMPONENTS	56
COATINGS ON CRYSTAL OPTICS	58
METALLIC COATINGS	60
· Front surface silver mirrors	60
· Front surface aluminum mirrors	62
· Special metallic coatings	64
SUBSTRATES	66
· High quality fused silica and CaF <sub>2</sub> components	66
· Substrate materials for UV, VIS and NIR/IR optics	68
· Transmission curves	69
· Standard substrates	70
· Prices	71
REGISTER	72

## FRONT SURFACE SILVER MIRRORS

### BROADBAND SILVER MIRRORS FOR THE VIS AND NIR




Figure 1: Reflectance spectra of silver mirrors with two types of protective coatings

**Optical properties:**

Standard:  $R \geq 94\%$  from 450nm to 600nm  
 $R \geq 95.5\%$  from 600nm to 1000nm  
 $R \geq 97\%$  for  $\lambda > 1\mu\text{m}$


Thin layer:  $R \geq 94\%$  from 450nm to 600nm  
 $R \geq 97\%$  from 600nm to 1000nm  
 $R \geq 98\%$  for  $\lambda > 1\mu\text{m}$

#### Special features:

- Highest reflectance of all metals in the VIS and NIR
- Sputtered silver mirrors show extremely low straylight losses ( $TS < 10^{-4}$  at 633nm)
- Silver must be protected because it is chemically unstable
- Sputtered protective layers guarantee very stable optical parameters and a lifetime of more than 4 years in normal atmosphere
- The thickness of the protective layer can be used to optimize the reflectance of the mirrors for different wavelength ranges (see example in figure 2).
- Optimized silver mirrors are ideal for use with femtosecond lasers (see pages 40–41)

#### Comparison of different types of protective layers


Property	Standard coating	Thin layer coating
Reflectance	Up to 2% lower than that of unprotected silver	Close to that of unprotected silver
Cleanability	o. k., tested according to MIL-M-13508C § 4.4.5	Cleaning is possible but the mechanical stability is worse than that of the standard coating





**Figure 2:** Reflectance spectrum of a silver mirror optimized for high reflectance between 1500 and 2000nm

## SILVER MIRRORS WITH ENHANCED REFLECTIVITY

The reflectivity of silver mirrors can be enhanced for selected wavelengths or wavelengths regions by a dielectric overcoat. Figures 3 and 4 show examples for such silver mirrors with enhanced reflectivity. Such mirrors combine very high reflectances at the wavelengths of interest with a relatively high reflectance throughout the VIS which makes them ideal for use with a pilot laser.



**Figure 3:** Silver mirror with enhanced reflectivity at 425 and 850nm ( $R > 99.5\%$ )




**Figure 4:** Reflectance spectra of a silver based scanning mirror for laser diodes in the NIR: HR ( $22^\circ - 58^\circ, 805 - 940\text{nm}$ )  $> 99.3\%$  combined with  $R > 50\%$  between 630 and 670nm

For more information on dielectrically enhanced silver mirrors see pages 12–15 and 41.

## FRONT SURFACE ALUMINUM MIRRORS

### BROADBAND MIRROR FOR THE UV, VIS AND NIR




**Figure 1:** Reflectance spectra of bare aluminum and of a standard protected aluminum mirror

#### Optical properties:

Bare aluminum:  $R > 80\%$  at 193 nm  
 $R = 92\%$  at 248 nm  
 $R > 85\%$  from 200 nm to 950 nm ( $R > 90\%$  from 230 nm to 600 nm)  
 $R > 90\%$  for  $\lambda > 1 \mu\text{m}$

Standard mirror:  $R = 82 \dots 92\%$  from 240 nm to 550 nm  
 $R = 85 \dots 92\%$  from 550 nm to 950 nm  
 $R > 92\%$  for  $\lambda > 1 \mu\text{m}$

The position of the minima and maxima of the reflectance depends on the design of the protective layer system and on the angle of incidence (AOI). Please specify AOI and the wavelengths of interest to optimize the reflectance as far as possible. The figure 2 shows the reflectance spectra of a mirror which was optimized for high reflectivity at 266 nm, 400 nm and 800 nm at AOI=45°.



**Figure 2:** Reflectance spectra of an aluminum mirror optimized for  $R > 85\%$  at 266 nm, 400 nm and 800 nm (AOI=45°, unpolarized light)

## FRONT SURFACE ALUMINUM MIRRORS

### OPTIMIZED ALUMINUM MIRRORS FOR THE DUV AND VUV




Figure 3: UV optimized Al: reflectance spectra of Al mirrors optimized for 157nm and 193nm (0°)

**Optical properties:** Special coating design depending on the wavelengths of interest

Optimized for 157nm:  $R=74 \dots 78\%$  for 157 nm ( $R>70\%$  from 150–200nm)

Optimized for 193nm:  $R=75 \dots 80\%$  for 193 nm

Optimized for 248nm:  $R>90\%$  for 248 nm

#### Special features:

- High reflectivity in the wavelength range specified
- Extremely low straylight losses ( $TS < 10^{-4}$  at 633 nm,  $TS < 10^{-3}$  at 248 nm,  $TS \approx 5 \times 10^{-3}$  at 193 nm)
- Standard mirrors can be cleaned using ethanol or acetone and are resistant to moderate abrasion (tested according to MIL-M-48497A § 4.5.4.2 and § 4.5.3.3)
- VUV optimized mirrors should be treated with extreme care
- All mirrors are resistant to humidity (tested according to MIL-M-13508C § 4.4.7)
- Highly stable optical parameters because of sputtered  $\text{SiO}_2$  protective layer

### ENHANCED ALUMINUM MIRRORS FOR 157nm

- Reflectance at 157nm can be further improved by dielectric overcoatings (up to  $R>94\%$ ).
- Reflectance in the VIS:  $R=60 \dots 80\%$ . This can be used for a pilot laser.
- Especially mirrors with  $R=85 \dots 90\%$  can be used at a wider range of AOI than all dielectric mirrors of this reflectivity.




Figure 4: Reflectance spectra of two types of enhanced Al mirrors for 157nm: a)  $R>80\%$  for AOI = 0° ... 45°, b)  $R>94\%$

## GOLD MIRRORS FOR THE NIR




Figure 1: Reflectance spectra of protected and unprotected gold mirrors

**Optical properties:** bare gold:  $R > 97\%$  from 700nm to  $1\ \mu\text{m}$   
 $R > 98\%$  for  $\lambda > 1\ \mu\text{m}$

Protected:  $R > 96\%$  from 700nm to  $2\ \mu\text{m}$   
 $R > 98\%$  for  $\lambda > 2\ \mu\text{m}$

**Special features:**

- High reflectance in the NIR
- Extremely low straylight losses ( $TS < 10^{-4}$  at 633 nm)
- Gold mirrors are chemically stable and can thus be used without protective layer
- Unprotected gold is soft and is easily scratched if it is touched
- Protected mirrors can be cleaned (tested according to MIL-M-13508C § 4.4.5)

## BROADBAND NEUTRAL DENSITY FILTERS FOR THE NIR/IR




Figure 2: Transmittance spectra of broadband neutral density filters with different transmittance values


- Broad wavelength range of nearly constant transmission ( $\sim 1\ \mu\text{m}$  to  $10\ \mu\text{m}$ )
- Substrate:  $\text{BaF}_2$
- Other transmittance values on request

## OTHER METALLIC COATINGS

- Chromium with extremely low pinhole density
  - Coatings on mask blanks and silicon wafers
- Zirkonium for use as filter for EUV radiation
- Molybdenum, platinum, copper, nickel and nickel/chromium according to customer specification

## METALLIC COATINGS – AN OVERVIEW

### Reflectance of metallic coatings



**Figure 3:** Reflectance of several metals versus wavelength (taken from A. McLeod "Thin film optical filters", A. Hilger, London, 1985)

- Silver: highest reflectance in the VIS and NIR, chemically unstable, soft, can only be used with protective layers, see separate data sheets on pages 40–41 and 60–61
- Gold: similar reflectance as silver in the NIR, chemically stable, soft, cannot be cleaned without protective layers
- Aluminum: relatively high and constant reflectance in the VIS and NIR, highest reflectance in the UV; surface oxide layer absorbs in the deep UV, soft, see separate data sheet on pages 62–63
- Chromium, platinum and rhodium can be used without protective layer

### Protective layers:

- Ensure cleanability and chemical stability
- Influence on the reflectance of the metal
- Very thin sputtered layers can be used for chemical protection of the metal because of high atomic density of the layers. Such layers show minimal influence on the VIS and NIR reflectivity of the metal.
- Even silver can be protected efficiently (lifetime more than 4 years for front surface mirrors).
- Mechanical protection and cleanability can only be reached by relatively thick protective layer systems.
- Optimization of the protective layer system for the wavelength of interest is particularly necessary in the UV.